




ESPECTROS DE INFRARROJO (IR)

PROBLEMA-PREPARACIÓN DE MUESTRAS




3048	39	2993	36	2948	21	2933	16	2863	26	2822	43	2747	66
2692	60	1689	4	1604	10	1572	6	1443	39	1421	49	1412	49
1392	63	1384	57	1372	55	1362	47	1347	33	1316	77	1263	8
1186	6	1162	20	1132	66	1031	79	963	18	956	35	923	62
888	62	874	50	824	66	732	79	718	77	712	74	697	70

ÁCIDO RETINOICO



ESPECTROS DE INFRARROJO (IR)

PROBLEMA-PREPARACIÓN DE MUESTRAS




3614	5	1511	26	1300	50	1149	47	843	23
3484	81	1494	18	1255	5	1103	8	588	84
3069	70	1477	47	1219	21	1044	58	527	50
3032	55	1464	28	1208	4	1038	52	438	70
2922	64	1443	55	1189	53	986	77		
2860	74	1381	68	1165	6	930	84		
1887	86	1325	10	1156	44	925	84		



ESPECTROS DE INFRARROJO (IR)

PROBLEMA-PREPARACIÓN DE MUESTRAS




3389	6	1418	37	1263	46	1070	4	770	49
3294	6	1361	49	1192	60	1046	9	741	57
2928	53	1342	44	1160	29	1036	8	703	58
2916	50	1309	49	1143	38	1020	20	669	49
2902	38	1290	62	1117	12	990	23	833	50
2878	43	1281	62	1105	14	892	55	604	41
1477	60	1266	60	1095	10	877	57	476	60

LACTOSA



ESPECTROS DE INFRARROJO (IR)

PROBLEMA-PREPARACIÓN DE MUESTRAS



3093	42	2606	77	1631	77	1236	9	826	68
3047	43	2487	81	1500	11	1168	41	812	31
3023	52	1934	81	1475	8	1153	52	753	9
2963	64	1846	84	1391	52	1072	49	690	12
2847	72	1606	22	1373	39	1024	64	618	64
2837	72	1598	4	1315	72	999	77	535	47
2723	74	1567	79	1293	77	888	64	607	38



ESPECTROS DE INFRARROJO (IR)


PROBLEMA-PREPARACIÓN DE MUESTRAS





ESPECTROS DE INFRARROJO (IR)

PROBLEMA-PREPARACIÓN DE MUESTRAS



3006	62	2589	66	1436	66	1222	21	918	26
2983	62	2548	68	1420	47	1190	4	841	50
2891	62	1754	9	1372	53	1136	58	805	52
2872	62	1693	6	1308	10	1096	66	766	41
2834	62	1606	14	1295	28	1014	57	706	46
2701	66	1483	64	1272	64	971	70	667	62
2670	66	1459	27	1257	64	928	62	600	68